




# Simultaneous Beam Sampling and Aperture Shape Optimization for Station Parameter Optimized Radiation Therapy (SPORT)

**Masoud Zarepisheh**

Ruijiang Li, Yinyu Ye, and Lei Xing



# SPORT Bridges the Gap




R. Li & L. Xing, Med Phys. 38, 4912-19, 2011.

## - Introducing an optimization framework to implement SPORT

- 1- R. Li and L. Xing, Medical Physics **38**(9), 4912–4919 (2011).
- 2- L. Xing, M.H. Phillips, and C.G. Orton, Medical Physics **40**(2), 020601 (2013).
- 3- R. Li and L. Xing, Medical Physics **40**(5), 050701 (2013).
- 4- L. Xing and R. Li, J. Phys.: Conf. Ser. **489**(1), 012065 (2014).




# Algorithm Flowchart





# Convergence: Head & Neck


## Head & Neck





# Convergence: Prostate

## Prostate






# SPORT V.S. IMRT

## Head & Neck

SPORT ———  
IMRT - - - -



# SPORT V.S. IMRT




## Head & Neck

Organ	Specific interest	SPORT	IMRT	Rel. diff.
PTV72	D (5%)	81.8	<b>81.6</b>	0.2%
	D (99%)	<b>69.4</b>	66.6	4.1%
PTV 54	D (5%)	<b>75</b>	75.9	1.2%
	D (99%)	<b>35.9</b>	32.2	10.8%
Cord	Max	<b>21.2</b>	30.9	37.2%
Brainstem	Max	<b>34.9</b>	39.4	12.1%
Mandible	Mean	31.9	<b>31.7</b>	0.6%
L-Parotid	Mean	13	13	0%
R-Parotid	Mean	<b>9.7</b>	11.4	16.1%



SPORT ———  
IMRT - - - -






# Optimized Angles V.S. Uniform

Head & Neck

SPORT-Optimized angles ———  
Uniform beam sampling - - -



# Optimized Angles V.S. Uniform



## Head & Neck


Organ	Specific interest	SPORT	Uniform Angles	Rel. diff.
PTV72	D (5%)	<b>81.8</b>	82	0.2%
	D (99%)	<b>69.4</b>	66.8	3.8%
PTV 54	D (5%)	<b>75</b>	<b>74.9</b>	<b>0.1%</b>
	D (99%)	<b>35.9</b>	34.7	3.3%
Cord	Max	<b>21.2</b>	32.6	42%
Brainstem	Max	<b>34.9</b>	37.5	7.1%
Mandible	Mean	<b>31.9</b>	33.3	4.3%
L-Parotid	Mean	<b>13</b>	15.3	16%
R-Parotid	Mean	<b>9.7</b>	14	36.3%



# Optimized Angles V.S. Uniform

Prostate

SPORT-Optimized angles —  
Uniform beam sampling - - -



# Summary & Future Research



- SPORT optimization framework with integration of three complementary algorithms
- Automatically determines the number of the stations & optimizes simultaneously the involved station parameters
- The optimal stations (control points) can be efficiently delivered with the new digital Linacs (.xml file & TrueBeam developer mode)
- SPORT with collimator modulation, non-coplanar and/or even no-isocentric beams

# Acknowledgments



- **Lei Xing**
  - Department of Radiation Oncology
- **Yinyu Ye**
  - Department of Management Science & Engineering
- **Ruijiang Li**
  - Department of Radiation Oncology